OTOMER

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EN _Wax plugs are due to the excessive secretion of ce ruminous glands. Endoscopic analysis shows that the wax appears to be a yellowish or brown mass partial ly or totally obstructing the external auditory canal. If wax plugs are not treated, they can even trigger the on set of serious diseases such as deafness, intense pruri tus, conductive hearing loss, etc.

OTOMER
Isotonic seawater solution specifically conceived to wash and clean adults' and children's ears to pre vent the formation of wax plugs and facilitate their extraction. The regular use of **Otomer** helps keep the external auditory canal clean and unobstructed, thus preventing the use of cotton swabs and any risk for the tympanum (eardrum).
Directions: remove the protective cap from the spray bottle and fit the nozzle on the valve. During distribu tion, the cylinder should be held in the vertical position to enable it to be completely emptied. Leaning head sli ghtly to one side, insert the nozzle into the external hea ring duct. Press the valve. Wait a few seconds then gent ly dry. Wash with hot water and dry the nozzle after use.
Frequency of use: use as directed, several times a day, if necessary.
WARNINGS: do not use in case of otitis, tympanum perforation, when foreign objects are obstructing the ear, or directly after ear surgery. No side effects or drug interactions have been reported. If product gets in the eyes, rinse thoroughly with water. Use under adult supervision.
Warning: pressurised container: may burst if heated. Keep away from heat, hot surfaces, sparks, open fla mes and other ignition sources. No smoking. Do not pierce or burn, even after use. Protect from sunlight. Do not expose to temperatures exceeding 50°C/122°F. Keep out of reach of children. Do not use other than as specified.

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EN_Ear hygiene spray

Medical Device

Read the instructions for use carefully, as they contain important information for the patient about the application of the medical device. Keep the instructions for further reference, if needed.

OTOMER Oil Spray is an effective, safe, simple product to soften and naturally remove ear wax and prevent water buildup in the ear while swimming (so-called "swimmer's ear"). Its use reduces the need for surgical removal of earwax by a doctor. When using OTOMER Oil Spray, the wax gradually dissolves and is naturally removed from the ear. Earwax is a secretion from glands located in the ear. Its excessive secretion can lead to the formation of deposits (plugs) of wax and as a consequence ear pain, discomfort and, in special cases, deterioration of hearing, including hearing loss.

Lingering water in the ear promotes bacterial growth and can lead to ear infections. OTOMER Oil Spray creates a protective barrier on the skin inside the ear, which effectively prevents water from accumulating inside the ear. The medical device OTOMER Oil Spray can be applied while standing, sitting or lying and disperses well into the ear canal.

WARNINGS AND PRECAUTIONS:

Pregnant women should seek medical advice before using the product. Always consult your doctor if you have frequent and unexplained outer ear disorders.

Do not use the medical device OTOMER Oil Spray:

in case of previously known hypersensitivity to any of the product components, if an eardrum perforation is found or suspected, in case of presence or suspicion of a foreign body in the ear, in case of ear leakage, without consulting a doctor in children under 3 months of age, without consulting a doctor in case of an earache accompanied by an elevated body temperature, in case of hearing disorders, in other body orifices, e.g. eye, nose.

Directions for use:

Medical device intended for administration in the ear.

1. Before first use, place applicator in a horizontal position and press dispenser 3-4 times to fill with liquid.
2. Gently pull the upper earlobe up and toward the back of the head.
3. Holding the bottle upright, insert the applicator into the external ear canal and press firmly on the top of the dispenser.
4. After application, massage the front of the ear canal from the outside.
5. If necessary, repeat steps 2 and 3 for the other ear.
6. After use, rinse dispensing tip in warm water and dry.
7. Turn applicator to vertical position to

prevent accidental dispensing of fluid.

Dosage:

Children over 3 months to 12 years of age:

- to remove wax buildup: 1 dose 2-3 times daily into each ear until wax is dissolved;
- for the prevention of earwax accumulation: 1 dose 1 time per week;
- as a preventive measure before bathing: 1 dose in each ear.

Use under adult supervision. Do not insert the applicator too deeply into the ear in order not to damage the auditory organ.

Children over 12 years of age and adults:

- to remove wax buildup: 1-2 doses 3 times daily into each ear until wax is dissolved;
- for the prevention of earwax accumulation: 1 dose 1 - 2 times per week;
- as a preventive measure before bathing: 1 dose in each ear.

If you have any doubt or feel that the effect of the medical device is too strong or too weak, contact your doctor or pharmacist.

Do not use higher doses than recommended without consulting your doctor.

Skipping a dose of the product has no significant effect on the course of therapy. Continue therapy at the current dosage. Do not use a double dose to make up for a missed dose.

Possible side effects

If symptoms of allergic reactions or other adverse reactions occur, discontinue use and consult a doctor if necessary. All adverse reactions should be reported to the manufacturer and to the distributor.

How to store: store in closed, original packaging at a temperature of 15-25° C. Keep out of reach and sight of children. Once opened, the product should be used within 6 months. Do not use medical device OTOMER Oil Spray after the expiration date stated on the package.

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ISTITUTO GANASSINI S.p.A.